

VOICE SWITCH

Version 1.1

An accessible voice activated switch using a HUZZAH32 – ESP32 Feather Board and an Amazon Echo.

Milador

Table of Contents

| Components | 2 |
|------------------|---|
| Hardware Setup | |
| Enclosure design | |
| Assembly | |
| Software Setup. | |

Github Repository: https://github.com/milador/VoiceSwitch

Components

- 1. HUZZAH32 ESP32 Feather Board (1 x)
- 2. ADAFRUIT NON-LATCHING MINI RELAY (1 x)
- 3. FEATHER HEADER KIT FML (1 x)
- 4. (OPTIONAL) BATTERY LITHIUM 3.7V 350MAH (1 x)
- 5. CONN JACK MONO 3.5MM PNL MNT (1 x)
- 6. HOOK-UP WIRE (1 x)
- 7. HEX NUT 0.197" NYLON M2.5 (2 x)
- 8. MACH SCREW PAN SLOTTED M2.5X0.45 (2 x)
- 9. (OPTIONAL) ECHO DOT OR OTHER ALEXA ENABLED DEVICE (1 x)

Hardware Setup


Enclosure design


The enclosure/housing files in STL format can be downloaded from GitHub repository under Hardware directory.

https://github.com/milador/VoiceSwitch/tree/master/Hardware/Enclosure


Please print the following STL files:


- 1 x VoiceSwitch_Top_Top.stl
- 1 x VoiceSwitch_Bottom.stl
- 1 x VoiceSwitch_Battery_Holder.stl (Optional)







Step 1: Position and solder the 12-position male header on HUZZAH32 board x 1.







Step 2: Position and solder the 16-position male header on HUZZAH32 board x 1.






Step 3: Position and solder the 3 pin Terminal Block of the Relay Feather Wing.







Step 4: Solder a wire from "Signal" pin to "A1" pin of the Relay Feather Wing.







Step 5: Position and solder the 14-position female header on the bottom of Relay Feather Wing.







Step 6: Position and solder the 16-position female header on the bottom of Relay Feather Wing.






Step 7: Solder a wire from "Sleeve" of 3.5mm Jack.






Step 8: Solder a wire from "Tip" of 3.5mm Jack.





Step 9: Solder the wire from "Tip" of 3.5mm Jack to "NO" of Terminal Block and "Sleeve" of 3.5mm Jack to "COM" of Terminal Block and.







Step 10: Place the Battery Holder in bottom enclosure. (Optional)





Step 11: Position the assembled HUZZAH32 board in bottom enclosure.






Step 12: Insert the M2.5 screws under the bottom enclosure and push them through the drill holes x 2.






Step 13: Insert the M2.5 screw nuts and secure HUZZAH32 board by tightening the nuts x 2.





Step 14: Stack the assembled Relay Feather Wing over the HUZZAH32 board.





Step 15: Insert the top enclosure over the bottom enclosure and snap them together.

Software Setup

- 1. Setup an Adafruit IO account according to the following instructions:
 - 1.1. Create an Adafruit IO account
 - 1.2. Create an AIO KEY
 - 1.3. Copy your *Adafruit IO* key for later use in the code
 - 1.4. Create an Adafruit IO feed with name of your switch. I used feed name *Monkey* as an example.
 - 2. Setup your HUZZAH32 ESP32 Feather Board
 - 2.1. Download ESP32 library for Arduino IDE and install it according to the following instruction: https://learn.adafruit.com/adafruit-huzzah32-esp32-feather/using-with-arduino-ide
 - 2.2. Download *Adafruit IO Arduino* library and add it to *libraries* directory under *Arduino* installation path in your computer.
 - https://github.com/adafruit/Adafruit IO Arduino
 - 2.3. Download Adafruit Arduino library for MQTT support and add it to libraries directory under Arduino installation path in your computer.
 - https://github.com/adafruit/Adafruit MQTT Library
 - 2.4. Download Adafruit Arduino HTTP Client library and add it to libraries directory under Arduino installation path in your computer.
 - https://github.com/arduino-libraries/ArduinoHttpClient
- 3. Setup your Alexa compatible device or Alexa application (Android and iOS) on your smartphone.
- 4. Setup IFTTT account 3.1. Create an IFTTT account at https://ifttt.com.
- 5. Create *Momentary switch* recipe
 - 5.1. Open *Ifttt* application on your smartphone.
 - 5.2. Click on the *Create* button to create a recipe.
 - 5.3. Click on *If This Add* button to add the trigger block.
 - 5.4. Search and select *Amazon Alexa* as trigger service.
 - 5.5. Click on *Say a specific phrase* block to select the trigger phrase.
 - 5.6. Enter the method phrase to activate the switch. I used phrase *Monkey switch* as an example.
 - 5.7. Click on *Then That Add* button to add the action service block.
 - 5.8. Search and select *Adafruit* as action service.
 - 5.9. Select **Send data to Adafruit IO** block.
 - 5.10. Enter *Monkey* feed as the *Feed Name*.
 - 5.11. Enter *switch* as the *Data to save*.

- 5.12. Click on the *Continue* button.
- 5.13. Click on the *Continue* button again.
- 5.14. Review the applet logic and click on the *Finish* button.
- 6. Create Switch on recipe
 - 6.1. Open Ifttt application on your smartphone.
 - 6.2. Click on the *Create* button to create a recipe.
 - 6.3. Click on If This Add button to add the trigger block.
 - 6.4. Search and select *Amazon Alexa* as trigger service.
 - 6.5. Click on *Say a specific phrase* block to select the trigger phrase.
 - 6.6. Enter the method phrase to activate the switch. I used phrase *Monkey on* as an example.
 - 6.7. Click on *Then That Add* button to add the action service block.
 - 6.8. Search and select *Adafruit* as action service.
 - 6.9. Select Send data to Adafruit IO block.
 - 6.10. Enter *Monkey* feed as the *Feed Name*.
 - 6.11. Enter *switch* as the *Data to save*.
 - 6.12. Click on the *Continue* button.
 - 6.13. Click on the *Continue* button again.
 - 6.14. Review the applet logic and click on the *Finish* button.
- 7. Create **Switch off** recipe
 - 7.1. Open *Ifttt* application on your smartphone.
 - 7.2. Click on the *Create* button to create a recipe.
 - 7.3. Click on *If This Add* button to add the trigger block.
 - 7.4. Search and select *Amazon Alexa* as trigger service.
 - 7.5. Click on **Say a specific phrase** block to select the trigger phrase.
 - 7.6. Enter the method phrase to activate the switch. I used phrase *Monkey off* as an example.
 - 7.7. Click on *Then That Add* button to add the action service block.
 - 7.8. Search and select *Adafruit* as action service.
 - 7.9. Select **Send data to Adafruit IO** block.
 - 7.10. Enter *Monkey* feed as the *Feed Name*.
 - 7.11. Enter *switch* as the *Data to save*.
 - 7.12. Click on the *Continue* button.
 - 7.13. Click on the *Continue* button again.
 - 7.14. Review the applet logic and click on the *Finish* button.

- 8. Setup the VoiceSwitch.ino code
 - 8.1. Open *VoiceSwitch.ino* file using Arduino IDE.
 - 8.2. Make a copy of your *IO_USERNAME* and *IO_KEY* from Adafruit IO account and paste them in the *config.h* file.
 - 8.3. Replace *your_ssid* with your WiFi's SSID and *your_pass* with your WiFi's password in the *config.h* file.
 - 8.4. Verify and upload *VoiceSwitch.ino* code to your HUZZAH32 ESP32 Feather Board.
- 9. Make sure the feed data status and status of your device match. You can manually set it to on/off to match it.